# --- STUDY PARAMETERS ---

| Date://202 Time::to:   |  |  |  |  |  |  |
|--|--|--|--|--|--|--|
| Researcher(s):    Site/Lab:  |  |  |  |  |  |  |
| PARTICIPANT  |  |  |  |  |  |  |
| Subject #: ID: Visit: Consent ☐ Safety ☐   |  |  |  |  |  |  |
| Height:cm Weight:kg Arm length:cm Ear-neck-fingertip:cm                          |  |  |  |  |  |  |
| Nasion-inion:cm L-R pre-auricular:cm Circumference:cm                            |  |  |  |  |  |  |
| Ethnicity: Medications: Nicotine: Alcohol:                                       |  |  |  |  |  |  |
| TMS  |  |  |  |  |  |  |
| Make: Magstim, MagVenture, Mag&More, Dantec, NeuroNav: 10/20, MNI, MRI,          |  |  |  |  |  |  |
| Model: 200, BiStim, Rapid   MagPro R, PowerMag   Other:                          |  |  |  |  |  |  |
| Coil shape: Round   Figure 8   Double-Cone   H-coil Handle: Flat   Branding iron |  |  |  |  |  |  |
| Coil size: 50   70   90   100   110   Other: mm Diameter: Inner   Mean   Outer   |  |  |  |  |  |  |
| EMG  |  |  |  |  |  |  |
| Muscle1: FDI,TE,ADM,FDS,EDC,FCR/U,ECR/U,BR,BB,TB,DEL,PEC,OtherL R                |  |  |  |  |  |  |
| Muscle2: FDI,TE,ADM,FDS,EDC,FCR/U,ECR/U,BR,BB,TB,DEL,PEC,OtherL R                |  |  |  |  |  |  |
| Muscle3: FDI,TE,ADM,FDS,EDC,FCR/U,ECR/U,BR,BB,TB,DEL,PEC,OtherL R                |  |  |  |  |  |  |
| Muscle4: FDI,TE,ADM,FDS,EDC,FCR/U,ECR/U,BR,BB,TB,DEL,PEC,OtherL R                |  |  |  |  |  |  |
| PNS  |  |  |  |  |  |  |
| Stimulator Make: Digitimer   ADInstruments Model: D   BioAmp   Other             |  |  |  |  |  |  |
| Nerve1: Digit_,Med,Uln,Rad, Site: Finger,Wrist,Elbow, Distance:cm L   F          |  |  |  |  |  |  |
| Nerve2: Digit_,Med,Uln,Rad, Site: Finger,Wrist,Elbow, Distance:cm L   R          |  |  |  |  |  |  |
| Nerve3: Digit_,Med,Uln,Rad,Site: Finger,Wrist,Elbow, Distance:cm L   R           |  |  |  |  |  |  |
| Nerve4: Digit_,Med,Uln,Rad, Site: Finger,Wrist,Elbow, Distance:cm L   R          |  |  |  |  |  |  |

## --- THRESHOLDS ---

#### TMS THRESHOLD LOCATIONS

| Site1:              | Ref.: Vertex,Inion,      | _ <b>Hem:</b> L,mid,R  | Right:cm              | Forward:c         | m O        |  |
|---------------------|--------------------------|------------------------|-----------------------|-------------------|------------|--|
| Site2:              | Ref.: Vertex,Inion,      | _ <b>Hem:</b> L,mid,R  | Right:cm              | Forward:c         | m O        |  |
| Site3:              | Ref.: Vertex,Inion,      | _ <b>Hem:</b> L,mid,R  | Right:cm              | Forward:c         | m <b>O</b> |  |
| Site4:              | Ref.: Vertex,Inion,      | _ <b>Hem:</b> L,mid,R  | Right:cm              | Forward:c         | m O        |  |
| TMS THRESHOLDS      |                          |                        |                       |                   |            |  |
| <b>Method</b> : MEP | Ps,Twitch,Other          | Algorithr              | n: Frequency,Pl       | EST,Other         |            |  |
| Target: 5/10/_      | MEPs of 10/20/ t         | rials <b>Criterion</b> | <b>mV:</b> 0.05/0.10/ | 0.20/1.0/Other_   | mV         |  |
| Muscle1:            | _ Side: L,R State: Rest/ | _%MVC,Move             | <b>mV</b> : 0.05/1    | Γhreshold:        | _%MSO      |  |
| Muscle2:            | _ Side: L,R State: Rest/ | _%MVC,Move             | <b>mV</b> : 0.05/1    | Threshold:        | _%MSO      |  |
| Muscle3:            | _ Side: L,R State: Rest/ | _%MVC,Move             | <b>mV</b> : 0.05/1    | Threshold:        | _%MSO      |  |
| Muscle4:            | _ Side: L,R State: Rest/ | _%MVC,Move             | <b>mV</b> : 0.05/1    | Γhreshold:        | _%MSO      |  |
| NERVE THRESHOLDS    |                          |                        |                       |                   |            |  |
| Nerve1:             | Side: L,R; Duration:r    | ms <b>Sensory:</b>     | mA; <b>M-wave:</b> _  | mA; <b>Twitch</b> | :mA        |  |
| Nerve2:             | Side: L,R; Duration:r    | ms <b>Sensory:</b>     | mA; <b>M-wave:</b> _  | mA; <b>Twitch</b> | :mA        |  |
| Nerve3:             | Side: L,R; Duration:     | ms <b>Sensory:</b>     | mA; <b>M-wave:</b> _  | mA; <b>Twitch</b> | :mA        |  |
| Nerve4:             | Side: L,R; Duration:r    | ms <b>Sensory</b> :    | _mA; <b>M-wave:</b> _ | mA; <b>Twitch</b> | :mA        |  |
|                     |                          | EEG                    |                       |                   |            |  |
| Hardware:           | Sampling ra              | ate:kHz F              | iltering, Low: _      | Hz <b>High</b> :_ | Hz         |  |
| Electrodes 1:       | :/2:/                    | 3:/                    | <b>4</b> :/           | Referen           | ce:        |  |
| Nerve:              | Side: L,R; Duration:ms   | s Intensity:           | _mA; Frequenc         | y:Hz Reps         | s:         |  |
| NOTES:              |                          |                        |                       |                   |            |  |
|                     |                          |                        |                       |                   |            |  |

## --- EXPERIMENT ---

### TMS LOCATIONS

| Site1: Hem.: L   Mid   R  |            |  |  |  |
|---|------------|--|--|--|
| Site2: Hem.: L   Mid   R Right:cm Forward:cm MNI:x,y,z                | : <b>C</b> |  |  |  |
| Site3: Hem.: L   Mid   R Right:cm Forward:cm MNI:x,y,z                | C          |  |  |  |
| Site4: Hem.: L   Mid   R Right:cm Forward:cm MNI:x,y,z                | : <b>C</b> |  |  |  |
| TMS PULSES  |            |  |  |  |
| TMS train type: Single   Paired   Triple   Quadruple   rTMS:Hz,s      |            |  |  |  |
| TMS coil 1: TMS coil 2: TMS coil 3: (in order of pulses)              |            |  |  |  |
| TMS intensity1: %MSO %RMT   %AMT Shape: Mono   Bi                     |            |  |  |  |
| TMS intensity2: %MSO %RMT   %AMT Shape: Mono   Bi                     |            |  |  |  |
| TMS intensity3: %MSO %RMT   %AMT Shape: Mono   Bi                     |            |  |  |  |
| TMS intensity4: %MSO %RMT   %AMT Shape: Mono   Bi                     |            |  |  |  |
| Total N blocks: N trials per block: N pulses per trial: Total pulses: |            |  |  |  |
| NERVE PULSES  |            |  |  |  |
| Nerve1: Digit, Med, Uln, Rad, Other Side: L,R Duration:ms Intensity:  | _mA        |  |  |  |
| Nerve2: Digit, Med, Uln, Rad, Other Side: L,R Duration:ms Intensity:  | _mA        |  |  |  |
| Nerve3: Digit, Med, Uln, Rad, OtherSide: L,R Duration:ms Intensity:mA |            |  |  |  |
| Nerve4: Digit, Med, Uln, Rad, Other Side: L,R Duration:ms Intensity:  | _mA        |  |  |  |
| Total N blocks: N trials per block: N pulses per trial: Total pulses: |            |  |  |  |
| NOTES   |            |  |  |  |
|   |            |  |  |  |
|   |            |  |  |  |
|   |            |  |  |  |
|   |            |  |  |  |